BISCO SPOTLIGHT

Marraskuu 2024

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1. Apply two separate coats of All-Bond Universal®, scrubbing the preparation with a microbrush for 10-15 seconds per coat. Do not light cure between coats.

2. Evaporate solvent by thoroughly air-drying with an air syringe for at least 10 seconds; there should be no visible movement of the material. The surface should have a uniform glossy appearance; otherwise, repeat Step 1-2. Light cure for 10 seconds.

3. Continue with placement of the restorative material according to the manufacturer's instructions.

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1. M.Sc.Dt. Emre ÖZEL, Dr. Haktan YURDAGÜVEN, Yrd.Doç.Dr. Esra CAN SAY, Prof.Dr. Sesin KOCAGÖZ, Evaluation of the Antibacterial Activity of Disinfectant Solutions with Phosphoric Acids Against Streptococcus Mutans. Journal of Hacettepe Faculty of Dentistry, Volume: 29, Issue 4, Page: 8-14, 2005
2. M. TURKUN1, Z. ERGUCU, L.S. TURKUN, E.U. CELIK, and M. ATES, Is Phosphoric Acid Sufficiently Antibacterial? J Dent Res 85 (Spec Iss B):abstract number 1605, 2006 (www.dental-research.org).

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1. Under rubber dam isolation, complete cavity preparation.

2. Achieve hemostasis. Leave dentin visibly moist.

3. Apply TheraCal LC® directly to the exposed pulp. Layer is not to exceed 1 mm in depth. Cover all the exposed areas and extend TheraCal LC at least 1mm onto sound dentin surrounding the exposure.

4. Light cure for 20 seconds. Place desired adhesive, base, and/or restoration following manufacturer's directions. Continue restoring tooth.

1. Isolate the tooth and perform conventional cavity preparation. Remove all infected carious tooth structure. Leave dentin visibly moist.

2. Apply TheraCal LC directly to the cavity floor of the preparation. Layer is not to exceed 1 mm in depth. Manipulate into a smooth surface covering all deep dentin areas.

3. Light cure for 20 seconds. Place desired adhesive, base, and/or restoration following manufacturer's directions. Continue restoring tooth.

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TAI KATSO
THERAFAMILY.COM

PAKKAUSKOOT

tuotenumero

37124 THERACAL LC, 1 g ruisku
37125 THERACAL LC, 4 x 1 g ruisku
36342 THERACAL LC -vientikärki 22 G, 50 kpl

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H-3301P
H-33014P
X-80621N

* BISCO has, on file, the calcium release data for TheraCal LC

1. Gandolfi MG, Siboni F, Prati C. Chemical-physical properties of TheraCal, a novel light-curable MTA-like material for pulp capping. International Endodontic Journal. 2012 Jun;45(6):571-9.

2. ADA definitions for direct and indirect pulp capping at: www.ada.org/en/publications/cdt/glossary-of-dental-clinical-and-administrative-terms

3. T. Okabe, M. Sakamoto, H. Takeuchi, K. Matsushima. Effects of pH on Mineralization Ability of Human Dental Pulp Cells. Journal of Endodontics. Volume 32, Number 3, March 2006.

4. Sangwan P; Sangwan A; Duhan J; Rohilla A. Tertiary dentinogenesis with calcium hydroxide: a review of proposed mechanisms. Int Endod J. 2013; 46(1):3-19

5. Selcuk SAVAS, Murat S. BOTSALI, Ebru KUCUKYILMAZ, Tugrul SARI. Evaluation of temperature changes in the pulp chamber during polymerization of light-cured pulp-capping materials by using a VALO LED light curing unit at different curing distances. Dent Mater J. 2014;33(6):764-9.

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						
BISCO Cements	TheraCem® Dual-Cured, Self-Adhesive Calcium & Fluoride Releasing Cement	Duo-Link Universal™ Dual-Cured Resin Luting Cement	Choice™ 2 Light-Cured Veneer Cement	BisCem® Dual-Cured, Self-Adhesive Resin Cement	eCement® LC Light-Cured Veneer Cement for Lithium Disilicate	eCement® DC Dual-Cured Resin Luting Cement for Lithium Disilicate
	Buy Now	Buy Now	Buy Now	Buy Now	Buy Now	Buy Now
Curing Mode	Dual Cured	Dual Cured	Light Cured	Dual Cured	Light Cured	Dual Cured
Requires Adhesive		●	●		●	●
Substrates	zirconia, glass ceramics, lithium disilicate, metal, alumina, composite	zirconia, glass ceramics, lithium disilicate, metal, alumina, composite	zirconia, glass ceramics, lithium disilicate, composite	zirconia, glass ceramics, lithium disilicate, metal, alumina, composite	lithium disilicate	lithium disilicate
Inlays/Onlays	●	●		●		●
Crowns/Bridges	●	●		●		●
Endodontic Posts	●	●		●		●
Maryland Bridges		●				●
Veneers			●		●	
Primer Required		●	●	●	●	●
Silane recommended for Glass Ceramics/Porcelain or Lithium Disilicate/Silicate	●	●	●	●	●	●
Best Uses	A self-adhesive cement, ideal for retentive preparations when "bonding" the restoration is NOT indicated. Releases both calcium and fluoride. Contains MDP.	Ideal for less retentive preparations when "bonding" the restoration IS indicated. Ideal for most restorations including CAD/CAM restorations.	Ideal cement for veneer restorations. Choice 2 is specifically formulated for high esthetics with color stability.	Translucent, dual-cured, self-adhesive resin cement. Indicated when an esthetic self-adhesive cement is preferred.	Designed for lithium disilicate veneer restorations (including CAD/CAM veneers).	Designed for lithium disilicate restorations (including CAD/CAM restorations).

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LUE LISÄÄ
THERACEM CA:STA

TAI KATSO
THERAFAMILY.COM

1. Gleave CM, Chen L, Suh BI. Calcium & fluoride recharge of resin cements. Dent Mater. 2016 (325):e26.
2. New Self-adhesive Resin Cement With Alkaline pH. Chen L, Gleave C, Suh B, J Dent Res96(A):#286, 2017

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

Itse levittyvä (Core-Flo DC Lite)

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

1. Isolate tooth and prepare accordingly. Rinse and dry.

2. Dispense 1-2 drops of **Universal Primer™ A & B** into a clean well. Apply two separate coats of **Universal Primer**, scrubbing the preparation with a microbrush for 10-15 seconds per coat.

3. Evaporate excess solvent by thoroughly air-drying with an air syringe for at least 10 seconds; there should be no visible movement of the adhesive. The surface should have a uniform glossy appearance; otherwise, apply an additional coat of **Universal Primer A & B** and air dry.

4. Place the end of the intra-oral tip directly onto the preparation, and express **Core-Flo™ DC** or **Core-Flo™ DC Lite**.

5. Light cure the **Core-Flo DC** or **Core-Flo DC Lite** for 40 seconds and continue with final preparation.

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

INTRAORAL REPAIRS OF PORCELAIN (PFM, OR PFZR)*

1. Isolate the area to be repaired. Remove the glaze and bevel (45 degrees) the porcelain around the area to be repaired. Sandblast or abrade with a coarse diamond bur. Rinse with water and air dry.

2. Place **Barrier Gel** on the gingival tissue and porcelain to protect areas where etching is not desired.

3. Apply 9.5% **Porcelain Etchant** to the dry porcelain surface for 90 seconds. Suction the etchant with high volume evacuation, then rinse with water and air dry. The etched porcelain should appear dull and frosty.

4. Apply 1 coat of silane (**Porcelain Primer**) to the etched porcelain and allow to dwell for 30 seconds. Dry with (warm) air syringe.

5. Apply 1 coat of **Z-Prime™ Plus** to the exposed metal/zirconia/alumina and dry with an air syringe for 3-5 seconds.

6. If metal masking is required, use **Dual-Cured Opaquer Base and Catalyst**, one drop of each, onto a mixing pad and mix with a brush tip. Apply a thin coat of the mixture only to the metal surface. Allow to self cure, or light cure for 5 seconds.

7. Apply a thin layer of **Porcelain Bonding Resin** to the repair site. Apply a hybrid composite (**Aelite™ All-Purpose Body**) to replace fractured porcelain and light cure, finish/polish.

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B-22300K

* Refer to instructions for complete details.

CONVENTIONAL CEMENTATION OF CROWNS (RETENTIVE PREPS)*

Using TheraCem

1. Sandblast internal surface of restoration (unless restoration has previously been sandblasted).

2. Remove temporary restoration and clean prep with pumice and water slurry. Rinse thoroughly and dry.

3. Try-in (to check occlusion and fit). After try-in, thoroughly rinse the restoration with water spray and dry.

4. Outside of the mouth, cover all bonded surfaces of the restoration with a layer of ZirClean®. Allow 20 seconds for the cleaning action of ZirClean to take effect, then thoroughly rinse with water spray and dry.

5. Dispense TheraCem® into the restoration. Seat the restoration with gentle passive pressure and remove any excess cement.

PRO TIP: To aid in the removal of excess cement, initially light-cure the margins for 2-3 seconds.

6. Light cure for 20-30 seconds or allow to self cure.

* Refer to instructions for complete details.

BONDED CEMENTATION OF INLAY/ONLAYS, BRIDGES, CROWNS (SHORT/TAPERED PREPS)*

Using Duo-Link Universal


OR

1a. Porcelain/Lithium Disilicate: Apply a thin coat of silane (**Porcelain Primer**) to the internal surface. Wait for 30 seconds, or dry with (warm) air.

1b. If restoration is metal/zirconia/indirect composite sandblast internal surface of restoration, unless restoration has previously been sandblasted.

2. Remove temporary restoration and clean prep with pumice and water slurry. Rinse thoroughly and dry.

3. Try-in (to check occlusion and fit). After try-in, thoroughly rinse the restoration with water spray and dry.

4. Outside of the mouth, Cover all bonded surfaces of the restoration with a layer of **ZirClean**®. Allow 20 seconds for the cleaning action of **ZirClean** to take affect, then thoroughly rinse with water spray and dry.

5. If restoration is metal/zirconia/indirect composite, apply one coat of **Z-Prime™ Plus** to the internal surface of the restoration and air dry for 3-5 seconds.

6. If prep includes enamel, selectively etch enamel with **Select HV® Etch w/BAC** for 15 seconds, suction then rinse thoroughly.

7. Dispense **Universal Primer™** and mix (or **All-Bond Universal®**) in a mixing well. Apply 2 separate coats, agitating each coat for 10-15 seconds.

8. Gently air dry until there is no visible movement of the adhesive. Then thoroughly air dry with greater air pressure. The surface should appear shiny, otherwise repeat step 7. Light cure for 10 seconds.**

9. Using a dual-cured resin cement (**Duo-Link Universal™**), fill the internal surface of the restoration and/or the prep with the cement. Seat with gentle, passive pressure and remove excess cement while holding restoration in place. Allow the cement to self-cure, or light-cure each surface of the tooth for 40 seconds.

* Refer to instructions for complete details.

** Universal Primer does not require light cure.

VENEER CEMENTATION*

Using Choice™ 2

OR

1. Prepare the internal surface of the restoration as instructed by the laboratory:

a. Porcelain/Lithium Disilicate:
If needed, etch the veneer with hydrofluoric acid (4% Porcelain Etchant or 9.5% Porcelain Etchant) according to the manufacturer's instructions.

b. Apply a thin coat of silane (Porcelain Primer/Bis-Silane™) to the internal surface. Wait for 30 seconds, or dry with (warm) air.

2. Remove the temporary restorations and clean the preparations (Cavity Cleanser™ and pumice slurry). Rinse thoroughly.

3. Try in the restorations using the corresponding shade of water-soluble try-in paste (Choice™ 2 Try-In Paste). Remove the veneers and either clean with etchant (Uni-Etch® w/BAC or Select HV® Etch w/BAC) and rinse thoroughly, or ultrasonicate in water or alcohol for 2-3 minutes. Dry the restorations.

4. Isolate. Etch preparations (Uni-Etch w/BAC or Select HV Etch w/BAC) for 15 seconds, rinse thoroughly.

5. Dispense adhesive (All-Bond Universal®) in a mixing well. Apply 2 separate coats, agitating each coat for 10-15 seconds.

6. Gently air dry until there is no visible movement of the adhesive. Then thoroughly air dry with greater air pressure. The surface should appear shiny, otherwise apply additional coats.

7. Light cure for 10 seconds.

8. Apply a thin layer of HEMA-free resin (Porcelain Bonding Resin) to the internal surface of the veneer. Do NOT light cure.


9. Using a light-cured resin cement (Choice 2), fill the internal surface of the restoration with the cement.


10. Seat with gentle, passive pressure and light cure for 2-3 seconds to tack the veneer into place. Remove excess cement, then light cure each veneer for 40 seconds.

* Refer to instructions for complete details.


VENEER CEMENTATION


CEMENTATION WITH A BONDING AGENT

CEMENTATION WITHOUT A BONDING AGENT


POST CEMENTATION WITH CORE BUILD-UP

Visit us at www.bisco.com
for a complete product
catalog! Rx only.

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